

Infection Control Division | Chemical Indicators

CG3

Three Line Labeler



Usage

Manual Labeler that allows a quick and easy labelling of sterilization packages with Chemdye® CD13 / CD23 / CD33 / CD43 / CD53 / CD83 Labels.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN, ISO 13485:2016.

Characteristics

Prints three lines of up to 12 characters each (Label size: 29 mm long x 22.5 mm wide)

Read instructions for use.

Accessory

Ink roller IRCG3.

Environmental conditions during manufacture

T = 15-30 °C, RH = 35-85 %.

Storage conditions

T = 10-30 °C, RH = 30-80 %.

Transport conditions

It is recommended to respect the storage conditions during transport.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

Warranty

The warranty coverage period is for one year from the purchase of the product.

Shelf life

Ink roller: 2 years.

Packing

1 labelling device and 1 ink roller per box.

Labelling

On product: product brand.

On labeler's box: product code and description, serial number, batch number, presentation, barcode and datamatrix code.

Infection Control Division | **Chemical Indicators**

CG3

Three Line Labeler



On ink roller: product code and description, batch number, presentation, expiration date and barcode.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical devices industries.

Other important information

Read product's instructions for use thoroughly before use.