

BS-H

Biological indicator holder

For Bionova® BT97 and BT94 biological indicators



Usage

Designed to hold BT97 and BT94 biological indicators during airborne and surface disinfection by Hydrogen Peroxide.

Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

Characteristics

High strength and durability aluminum re-usable holder.

Dimensions: 54 mm x 20 mm x 21 mm

Place the holder with the indicator in those areas that you consider a priori more inaccessible to the Hydrogen Peroxide.

Storage conditions

T = 10-30 °C, RH= 30-80 %, store in dark place.

Transport conditions

Carry in closed and strengthened boxes in order to avoid damages. The transport of this product does not represent a risk for health.

Packing

1 unit per box.

Packing information: product description, manufacturer information, storage conditions and presentation

Labelling

On product's packing: product code, batch number, manufacture date, bar code and datamatrix code.

Possible target markets

Healthcare, Food, Pharmaceutical and Medical industries.

Other important information

Read Bionova® BT97 and BT94 instructions for use thoroughly before use.

Precautions

Store the product within its' original package, away from disinfectants, corrosive or oxidant agents, gases and humidity.