

# IRCG3

## Ink Roller



### Usage

Spare Ink Roller is made to fit the CG3 – Three Line Labeler Chemdye® Terragene®, which allows a quick and easy labelling of sterilization package with Chemdye®CD13 / CD23 / CD33 / CD43 / CD53 / CD73 / CD203 / CD403 Labels.

### Applicable Regulation

N/A

### Product Classification

Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012.

### Characteristics

This spare roller contains black ink for a consistent label every time. The roller is designed to be used with CG3 – Three Line Labeler. This replacement ink roller is economical to keep plenty on hand. This long-lasting roller is easily released from the labeler and snaps into it quickly and simply, making previous roller change easy. Read CG3 directions for use.

•PERFORMANCE: 1 box x 12 rolls of labels (750 labels per roll).

**WARNING:** performance will depend on the pressure exerted on the labeler at the time of printing each label (which will alter the amount of ink that is printed on each label), as well as the number of digits per line used for printing

### Storage conditions

T = 10-30 °C, RH 30-80 %.

### Transportation conditions

Transport in closed and reinforced boxes in order to avoid damages. Transportation of this product does not represent a risk for health.

### Shelf-life

N/A

**Packing**

1 ink roller per package.

**Labelling**

On product's package: product code, description and batch number.

**Possible target markets**

Healthcare, Food, Pharmaceutical and Medical industries.

**Other important information (optional)**

Read CG3 directions for use thoroughly before use.

**Note:** When necessary, limits for each product's parameters should be specified.